

ENG INSTRUCTION MANUAL FOR DIGITAL BLOOD PRESSURE MONITOR LD
POL INSTRUKCJA OBSŁUGI CIŚNIENIOMIERZA CYFROWEGO LD
HUN HASZNÁLATI ÚTMUTATÓ LD KÉSZÜLÉK A VÉNYOMÁS ÉS PULZUSSZÁM DIGITÁLIS MÉRÉSERE
ROU MANUAL DE UTILIZARE A TENSIOMETRULUI DIGITAL LD
BGR ПРКЛОДВОБОДО ЗА ЕКСПЛУАТАЦИЯ НА ЦИФРОВ ТОНОМЕТР LD

fig.1 rys. 1 1. ábra desenul 1 рис. 1
 PARTS AND COMPONENTS DENUMIREA PĂRȚILOR SI COMPONENTELOR НАЗВАНИЯ ЧАСТИ И КОМПОНЕНТОВ
 PODSTAVOWE ČĘŚCI I KOMPONENTY НАЗВАНИЯ ЧАСТИ И КОМПОНЕНТОВ
 АЛКАТРАЗЕК МЕНЕВЕЗЕ




fig.2 rys. 2 2. ábra desenul 2 рис. 2
 BATTERY INSTALLATION WYMIAŁA BATERII INSTALAREA BATERIILOR ПОСТАВЉЕЊЕ НА ЕЛЕМЕНТИТЕ НА ЗАХРАНВАЊЕ




fig.3 rys. 3 3. ábra desenul 3 рис. 3
 CORRECT POSITION DURING MEASUREMENT PRZYJĘCIE POZYCYJI UMOŻLIWIĄcej POMIAR A HELYES TESTELHETZET A MÉRSÉLN POSZTÍTÁSA NELMÉN PÓZITA CORECTĂ ÎN TIMPUL PROCESULUI DE MÁSURARE ПРАВИЛНА ПОЗА ПРИ ИЗМЕРВАЊЕ




fig.4 rys. 4 4. ábra desenul 4 рис. 4
 CUFF PREPARATION ZAKŁADANIE MANKIETU A MANDZSETTA ELOKÉSZÍTÉSE




fig.5 rys. 5 5. ábra desenul 5 рис. 5
 MEASUREMENT PROCEDURE PROCEDURA POMIARU A MÉRSÉ MODSZERE





fig.6 rys. 6 6. ábra desenul 6 рис. 6
 MÉMORY FUNCTION FUNKCIJA PAMIĘCI FUNKCIJA, „ПАМЕТ“

**ENG**

DENOMINATION OF PARTS AND COMPONENTS
 (Fig. 1)
 ① Electronic block
 ② LCD display
 ③ Cuff socket
 ④ Cuff connector
 ⑤ Air hose
 ⑥ Cuff (applied part)
 GENERAL INFORMATION

This manual is intended to provide assistance to user in terms of safe and effective operation of arterial blood pressure and heart rate measuring device with value display led (version LD-521, LD-521A) (hereinafter referred to as the "DEVICE"). The Device should be used in accordance with the rules stated in this manual and should not be used for the purposes other than those described herein. It is important to read and understand the entire manual, especially the section "Recommendations for proper use".

INDICATIONS FOR USE

The device is intended for measurement of systolic and diastolic arterial blood pressure and determining pulse rate in patients at the age of 15 and up. The device is recommended for use by patients with unstable (inconstant) arterial blood pressure or with known arterial hypertension at home in addition to medical supervision.

PRINCIPLE OF OPERATION

Oscillometric method of arterial blood pressure and heart rate measurement is used in the device. The cuff is wrapped around the upper arm and pumped automatically. The sensitive element inside the device detects weak fluctuations of pressure in the cuff generated by expansion and contraction of the brachial artery in response to each heart beat.

Amplitude of the pressure waves is measured, converted into millimeters of mercury column and is displayed in the form of a digital value. The device has 90 memory cells to store the measurement results. Please note that the device may not ensure the specified measurement precision, if it is used at temperature or moisture other than those stated in the section "Technical characteristics" of this manual. We give notice of possible errors in measurement of arterial blood pressure in this device with pronounced arrhythmia. Please consult your physician as to how to measure arterial blood pressure in children.

The patient is an intended operator. But persons who suffer from arrhythmia, diabetes, cardiovascular problems or who have had a stroke should consult your doctor before using the device.

IMPORTANT SAFETY INSTRUCTIONS

- Do not use the device with the cuff fully in use with patient.
 - Do not use the device with other medical electrical (ME) equipment simultaneously.
 - Do not use the device in the area of the HF surgical equipment, MRI, or CT scanner exists, or in the oxygen rich environment.

- Do not use a mobile phone or any other devices that emit electromagnetic fields, near the device. This may result in incorrect operation of the device.
 - Never use any accessories or parts from other manufacturers. Using such accessories or parts could cause a hazardous situation for the user or damage to the device.
 - Do not modify this equipment without authorization of the manufacturer.

RECOMMENDATIONS ON CORRECT MEASUREMENT

1. For correct measurement you should know that THE BLOOD PRESSURE IS SUBJECT TO SHARP VARIATIONS EVEN WITHIN THE SHORT TIME INTERVALS. The blood pressure depends on many factors. It is usually lower in summer and higher in winter. The blood pressure varies together with the atmospheric pressure, depends on physical loads, emotional excitement, stress and dietary regime. Drugs, drinking alcohol and smoking produce significant effect on blood pressure. Even the very procedure of blood pressure measurement in a hospital sends the blood pressure high in many people, thus the blood pressure measured at home often differs from the values received in a hospital.

As the blood pressure tends to rise at low temperatures, make measurements at an indoor temperature (approximately 20°C). If this Device was stored in the room with low temperature, keep it at least 1 hour at indoor temperature before use, otherwise the measurement result may be incorrect. During the difference in readings for healthy people may be 30-50 mmHg of systolic pressure and to 10 mmHg of diastolic pressure. The dependence of the blood pressure on various factors is individual for each person. Accordingly, it is recommended to keep a special book with blood pressure measurement records. ONLY A CERTIFIED DOCTOR IS CAPABLE TO ANALYZE THE TENDENCY OF YOUR BLOOD PRESSURE VARIATIONS USING YOUR RECORDS.

2. People with cardiovascular and some other diseases requiring blood pressure monitoring should make measurements in the hours fixed by attending doctor. REMEMBER THAT THE DIAGNOSTIC AND ANY TREATMENT OF HYPERTENSION MAY BE CONDUCTED ONLY BY A CERTIFIED DOCTOR ON THE BASIS OF BLOOD PRESSURE VALUES OBTAINED BY THIS DOCTOR. DRUGS PRESCRIPTION SHOULD NOT BE CONTROLLED ONLY BY YOUR ATTENDING DOCTOR.

3. A person with cardiovascular disease, weak pulse and also in patients with the prominent distortions of cardiac rhythm it is not difficult to determine blood pressure accuracy. In SUCH CASES CONSULT A CERTIFIED DOCTOR ABOUT APPLICATION OF THE ELECTRONIC DEVICE.

The following persons should also consult their doctor before using the device: persons with abnormalities in the shoulder area that interfere with normal blood flow; persons with disorders of the peripheral circulatory system (diabetes mellitus, cirrhosis, atherosclerosis, etc.); persons who use implantable or wearable medical electrical equipment, such as a pacemaker, defibrillator, or electrocardiograph monitor; persons undergoing hemodialysis or treatment with anticoagulants, anti-platelet drugs, or steroids. THE PULSE RATE INDICATOR IS NOT INTENDED FOR PACEMAKERS' OPERATION MONITORING.

4. KEEP QUIET DURING A MEASUREMENT TO OBTAIN THE ACCURATE VALUES OF YOUR BLOOD PRESSURE WITH THE ELECTRONIC DEVICE. Measure your blood pressure in the calm and comfortable conditions at the indoor temperature. No eating an hour before measurement; no smoking, no taking tonic agents or alcohol 1.5 hours before measurement.

5. The accuracy of blood pressure measurement depends on whether the cuff matches the size of your arm. THE CUFF SHOULD NOT BE TOO SMALL OR TOO LARGE.

6. Wait 3 minutes between measurements for the blood to restore its circulation. However, the persons with prominent atherosclerosis due to considerable loss of vascular elasticity must decrease the wait time between measurements (10-15 minutes). This also refers to the patients suffering from diabetes. For more accurate determination of blood pressure, it is recommended to make a series of 3 consecutive measurements and to use the average value.

POWER SUPPLY OF THE DEVICE
 Installation of batteries (fig. 2)

1. Open the cover of the battery compartment and insert 4 batteries of type AAA as shown on the figure inside the compartment. Make sure the polarity is correct. Do not apply undue exertion when opening the cover of the battery compartment.

2. Close the battery compartment cover.
 • Replace all batteries when the indicator of batteries replacement is being constantly displayed or there is no indication on the display. The battery replacement indicator does not show the level of discharge.

The batteries that come with the device are intended for check of the device serviceability when purchased, their service life can be less than that for recommended batteries.

• When replacing the batteries replace all of them simultaneously. Do not use second-hand batteries.
 • If the device is not used for long time, remove the batteries from the device.
 • Do not leave used batteries in the device.

Use of the device with AC adapter
 The AC adapter is specified as a part of the blood pressure monitor. Use the AC adapter (model LD-N063, supplied as part of package for model LD-521A) with the following technical characteristics.

Operating voltage: 3V ± 5%; 500 mA
 Type of connector: micro-USB

The socket for the stabilized power supply is located on the right side of the device.

In order to use the device with a power source connect the power source connector to the device and insert the power source plug into the socket outlet and press the button . Upon completion of measurement switch the device off by pressing the button , remove the power source plug from the socket outlet and disconnect the power source connector from the device.

CORRECT POSITION DURING MEASUREMENT (fig. 3)

1. Sit at a table so that your arm rests on its surface during measurement. Make sure that the point of the cuff application is approximately at the same height as your heart, and that your forearm lies freely on the table and does not move.

2. You can measure pressure when lying on your back. Look at the ceiling, keep calm and do not move during measurement. Make sure that the measurement point on your arm is approximately at the same level as your heart.

CUFF PREPARATION (fig. 4)

1. Insert the cuff end for about 5 cm into a metal ring as shown in figure.

2. Apply the cuff to your upper arm so that the air tube is directed to your palm. If the measurement on your left arm is difficult, you may use your right. In this case remember that the readings may differ by 5-10 mmHg and even more.

3. Wrap the cuff around your upper arm so that the bottom of the cuff is approximately 2-3 cm above your elbow. The sign should be over your heart.

4. Fix the cuff so that it fits tightly on the arm, but make sure that it is not overtight. Too tight or too loose placement of the cuff may lead to inaccurate measurement.

5. On the front of the cuff the sign should point to the area . It means that the cuff is chosen correctly and fits to the size of your upper arm. If the sign points to the area marked the cuff is too small and the readings will be higher. If the sign points to the area marked the cuff is too large and the readings will be lower.

6. If the arm has a coniform, the cuff should be put with a spiral movement as shown in the figure.

7. If the rolled-up sleeve squeezes the arm interfering with free blood flow the device may give inaccurate figures not corresponding to your actual blood pressure.

MEASUREMENT PROCEDURE (fig. 5)

1. Insert the power connector into the cuff socket. Before measurement take a deep breath 3 to 5 times and relax. Do not move, do not talk and do not strain your arm during measurement.

2. Press the button .

3. All signals will appear for a short time on the display (Fig. 5.1), a short signal will sound and the air will start quickly pump into the cuff.

4. After reaching 120-130 mmHg the cuff will start pumping additionally, first slowly, then quickly. The pulse symbol (Fig. 5.2) will start flashing during pumping. The automatic inflation Measuring System (IMS) allows determination of systolic and diastolic blood pressure in the process of inflation. SINCE THE ARTERIAL BLOOD PRESSURE AND THE PULSE ARE MEASURED DURING INFLATION OF THE CUFF, TRY TO REMAIN IMMOBILE AND DO NOT MOVE YOUR ARM DURING MEASUREMENT AS WELL AS NOT TO STRAIN MUSCLES OF YOUR ARM.

5. After the end of measurement a signal will be heard, after which the device will release air from the cuff and the measurement result will appear on the display (Fig. 5.3). The flashing symbol on the display informs about irregular pulse. Display of arrhythmia indicator may be caused by body movement during measurement. Consult your physician, if you regularly observe the symbol .

6. The pressure in the cuff will drop rapidly. The automatic inflation Measuring System (IMS) allows determination of systolic and diastolic blood pressure in the process of deflation. SINCE THE ARTERIAL BLOOD PRESSURE AND THE PULSE ARE MEASURED DURING DEFLATION OF THE CUFF, TRY TO REMAIN IMMOBILE AND DO NOT MOVE YOUR ARM DURING MEASUREMENT AS WELL AS NOT TO STRAIN MUSCLES OF YOUR ARM.

7. The device will release air from the cuff and the measurement result will appear on the display (Fig. 5.4).

8. The measurement procedure is completed.

FORCED DECOMPRESSION OF THE CUFF

If during inflation of the cuff you need to quickly decompress the cuff, press the button .

The device will release air from the cuff and the measurement result will appear on the display (Fig. 5.5).

MEMORY FUNCTION (fig. 6)

1. When you press the button the device will display the last 9 measurements.

2. The device will display the last 9 measurements.

3. The device will display the last 9 measurements.

4. The device will display the last 9 measurements.

5. The device will display the last 9 measurements.

6. The device will display the last 9 measurements.

7. The device will display the last 9 measurements.

8. The device will display the last 9 measurements.

9. The device will display the last 9 measurements.

10. The device will display the last 9 measurements.

11. The device will display the last 9 measurements.

12. The device will display the last 9 measurements.

13. The device will display the last 9 measurements.

14. The device will display the last 9 measurements.

15. The device will display the last 9 measurements.

16. The device will display the last 9 measurements.

17. The device will display the last 9 measurements.

18. The device will display the last 9 measurements.

19. The device will display the last 9 measurements.

20. The device will display the last 9 measurements.

21. The device will display the last 9 measurements.

22. The device will display the last 9 measurements.

23. The device will display the last 9 measurements.

24. The device will display the last 9 measurements.

25. The device will display the last 9 measurements.

26. The device will display the last 9 measurements.

27. The device will display the last 9 measurements.

28. The device will display the last 9 measurements.

29. The device will display the last 9 measurements.

30. The device will display the last 9 measurements.

31. The device will display the last 9 measurements.

32. The device will display the last 9 measurements.

33. The device will display the last 9 measurements.

34. The device will display the last 9 measurements.

35. The device will display the last 9 measurements.

36. The device will display the last 9 measurements.

37. The device will display the last 9 measurements.

38. The device will display the last 9 measurements.

